

Avoid Common Missteps

Scan the QR code below



or visit jnvisionpro.com/calculators-tools
to **QUICKLY AND EASILY** find the right trial
lens powers for presbyopic patients

Keep these calculators handy on your device!

Save fitting calculator to your desktop or
mobile device home screen

1. JJV Data on File 2020. ACUVUE® PUPIL OPTIMIZED DESIGN TECHNOLOGY: JJVC Contact Lenses, Design Features, and Associated Benefits.

The calculators are provided to clinicians as a convenience, in conjunction with the labeling of ACUVUE® MULTIFOCAL Contact Lenses and are not a substitute for professional judgment.

Important Information: ACUVUE® Contact Lenses are indicated for vision correction. As with any contact lens, eye problems, including corneal ulcers, can develop. Some wearers may experience mild irritation, itching or discomfort. Lenses should not be prescribed if patients have any eye infection, or experience eye discomfort, excessive tearing, vision changes, redness or other eye problems. Consult the package insert for complete information. Complete information is also available from Johnson & Johnson Vision Care, Inc. by calling 1-800-843-2020, or by visiting www.jnvisionpro.com.

©Johnson & Johnson Vision Care, Inc. 2022 | PP2020MLT8266v2.0 | MB-00637

take **COMFORT** in
ACUVUE®

AOMF72206

Fit Success Tips

ACUVUE®
MULTIFOCAL
WITH PUPIL OPTIMIZED DESIGN

The only brand with 100% of parameters
optimized by both age and refraction
vs. <2% for the leading competitor*¹

* Compared to the leading competitor's design; technology optimized for both the parameters of refractive error and add power.

Fit in 3 Easy Steps

1 Perform fresh spherical refraction in ambient light

- **Avoid over-minusing** by performing a red/green balance
- **Apply vertex correction and determine spherical equivalent** (or, use the digital fitting calculator)
- Ensure cylinder $\leq 0.75D$

2 Determine ocular dominance using sensory method

- **Utilize Sensory Method** (+1.50 blur tolerance test)
- Remember, sensory dominance and sighting dominance may **differ up to 40% of the time¹**

3 Determine functional add

- **Prescribe minimum plus** for functional near vision
- **Avoid prescribing additional plus** in the ADD, as is occasionally done in progressive spectacle lenses

Find the Right Trial Lens

• Scan the QR code on the back, or to go to the simple online fitting calculator at jnvisionpro.com/calculators-tools

• If the fitting calculator is not used:

- Refer to the fit selection table below based on the ADD to determine initial trial lenses
- If your patient still has a specific near or distance need, refer to the enhancement tables below under "Problem Solving"

Finding the right lens for the right balance of distance and near vision

ADD	EYE	LENS SELECTION
+0.75D to +1.25D	Dominant Eye	LOW
	Non-dominant Eye	LOW
+1.50D to +1.75D	Dominant Eye	MID
	Non-dominant Eye	MID
+2.00D to +2.50D	Dominant Eye	MID
	Non-dominant Eye	HIGH

Problem Solving: Finding the right lens to adjust for enhanced distance vision or enhanced near vision

ENHANCED DISTANCE VISION	ENHANCED NEAR VISION
ACUVUE® SPHERICAL CONTACT LENS	LOW
LOW	LOW +
LOW	MID
MID	MID +
MID	MID
MID +	HIGH +

Binocular Over-Refraction to check distance power

To Improve Distance:

- For Low & Mid ADD, lower ADD in dominant eye
- For High ADD, lower ADD in non-dominant eye and add +0.25D

To Improve Near:

- ADD +0.25D to non-dominant eye

+ Add +0.25D to the distance power

10
min

Allow for 10 minutes of real-world exposure (outside of the exam room) before assessing visual performance.